
cobas[®] PCR sample kits

The first step to an accurate CT/NG result


A short guide for collecting & transporting samples for testing

Collecting female endocervical swab samples


1. CLEAN: Using one of the swabs (provided), remove excess mucus from the cervical os and surrounding mucosa. **Discard this swab after use.**


Note: Cleaning excess mucus from the cervical os is required to assure adequate target cells are obtained for processing. A large-tipped cleaning swab may be used (not provided).

2. COLLECT: To collect the specimen, insert the second swab provided into the endocervical canal. Gently rotate the swab 5 times in one direction in the endocervical canal. Do not over rotate. Carefully withdraw the swab, avoiding any contact with the vaginal mucosa.

3. ALIGN: Remove the cap from the **cobas[®]** PCR Media tube and lower the swab specimen into the tube until the visible dark line on the swab shaft is aligned with the tube rim. The tip of the swab should be just above the media surface near the hexagonal Roche logo.

4. BREAK: Carefully leverage the swab against the tube rim to break the swab shaft at the dark line; discard the top portion of the swab.

5. CAP: Tightly re-cap the **cobas[®]** PCR Media tube.

Handling precautions:

- **Do not** pre-wet collection swabs with the collection media before obtaining the endocervical specimen.
- Use care to avoid splashing of contents.

cobas[®] PCR Media tube: Shipping and storage

- After collection, you can ship your **cobas[®]** PCR Media vials at ambient temperature with the specimen collection swab inside.
- Endocervical swab specimens collected in **cobas[®]** PCR Media are stable at 2 - 30° C for up to 1 year.

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